

THRESHOLD RAMP



Ideal for sliding doors (two ramps shown)

The EZ-ACCESS® Modular Threshold ramp is a lightweight yet durable free-standing ramp designed for doorways, sliding glass doors and raised landings. The modular design, with its interchangeable components, easily accommodates thresholds from ¾" to 6" high. This seam-free ramp is 34" wide and ideal for use with wheelchairs, walkers, and cane users. Made of anodized aluminum and is intended for indoor or outdoor use.



Modular design features an adjustable top flap

SPECIFICATIONS	THRESH 1	THRESH 1.5	THRESH 2	THRESH 3	THRESH 4	THRESH 5	THRESH 6
Made of aluminum	•	•	•	•	•	•	•
Non-skid surface	•	•	•	•	•	•	•
Usable size (H x L x W)	1" x 6½" x 34"	1½" x 9½" x 34"	2" x 12" x 34"	3" x 17" x 34"	4" x 22¼" x 34"	5" x 27½" x 34"	6" x 32¾" x 34"
Weight (lbs.) *	3	5	5½	8	11	15	20
Weight capacity (lbs.)	600	600	600	600	600	600	600
For wheelchairs & scooters	•	•	•	•	•	•	•
Usage guidelines	¾" - 1¼" rise	1¼" - 1¾" rise	1¾" - 2¼" rise	2¾" - 3¼" rise	3¾" - 4¼" rise	4¾" - 5¼" rise	5¾" - 6¼" rise

* Weights are approximate and may vary.

Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.

