



Redesigned Top to Bottom

We are proud to introduce the SUITCASE® Signature Series® portable wheelchair ramp. Its single-fold design sets up quickly and folds in half lengthwise to be carried like a suitcase when not in use and features several innovative design enhancements. From the no pinch, non-protruding full-length live hinge and durable, ergonomically designed handles to the enhanced extruded tread for superior traction, this ramp was redesigned from the ground up to be simply the best in the industry. Top Lip Extension (TLE) compatible. Available in 2' – 8' lengths with an 800-pound weight capacity.



Features a no-pinch, non-protruding live hinge for safe, silent, and maintenance free operation




Ergonomically designed flexible, non-breakable handles offer convenient carrying and a comfortable grip



Bottom transition plates self-adjust to accommodate uneven terrain

OPTIONAL TOP LIP EXTENSION (TLE)

The TLE will increase the reach of the SUITCASE® Signature Series® portable wheelchair ramp by extending the lip from the standard 3" length to 9". The EZ-ACCESS TLE is unique because it can be added or removed as needed, eliminating the need for a large ramp inventory. Shipped as a pair.



With the use of a TLE, the rear bumper is easily cleared

SPECIFICATIONS	SUITCASE SS2	SUITCASE SS3	SUITCASE SS4	SUITCASE SS5	SUITCASE SS6	SUITCASE SS7	SUITCASE SS8	TLE
Made of aluminum	•	•	•	•	•	•	•	•
Non-skid surface	•	•	•	•	•	•	•	•
Usable size	2' x 30"	3' x 30"	4' x 30"	5' x 30"	6' x 30"	7' x 30"	8' x 30"	8½" x 13 ⁷ / ₈ " x 2 ³ / ₈ "
Folded size (L x W x H)	24" x 15" x 4"	36" x 15" x 4"	48" x 15" x 4"	60" x 15" x 4"	72" x 15" x 4"	84" x 15" x 4"	96" x 15" x 4"	n/a
Built-in carrying handle	•	•	•	•	•	•	•	n/a
Weight (lbs.) *	9	14	18	25	30	36	42	5 per set
Weight capacity (lbs.)	800	800	800	800	800	800	800	**500
For wheelchairs & scooters	•	•	•	•	•	•	•	•
Usage guidelines	Curbs, 1 step	Curbs, 1 step	Curbs, 1-2 steps	Curbs, 1-2 steps, some minivans	1-3 steps, some minivans & SUVs	1-3 steps, some minivans & SUVs	2-3 steps, some minivans & SUVs	n/a

* Weights are approximate and may vary. Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.

** Ramp weight capacity is reduced when used with the TLE.



Bariatric weight capacity = 800 lbs. or greater



UPS Shippable